

## 50 μL, Model 1705 RN SYR, Small Removable NDL, 22s ga, 2 in, point style 2

Specification Sheet (Part/REF # 80930)



Split/splitless injection consists of a heated chamber with a glass liner into which the sample is injected through the septum. A microsyringe is used to inject the sample through a rubber septum into a flash vaporizer chamber at the head of the column.

## **Product Specifications**

| Volume                | 50 μL                         |
|-----------------------|-------------------------------|
| Syringe Type          | Gastight                      |
| Syringe Series        | 1700                          |
| Termination           | Removable Needle (RN)         |
| Needle                | Small RN                      |
| Gauge                 | 22s gauge                     |
| Point Style           | 2                             |
| Needle Length         | 2 inch (51 mm)                |
| Fluid Path            | SST, borosilicate glass, PTFE |
| Autoclavable          | Yes                           |
| Gas Sterilizable      | Yes                           |
| Pressure              | 1000 psig (68.9 bar)          |
| Maximum Temperature   | 115°C / 239°F                 |
| Minimum Temperature   | 10°C / 50°F                   |
| Barrel Inner Diameter | 0.041 inches (1.03 mm)        |
|                       |                               |



Visit <u>www.hamiltoncompany.com/support</u> to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055



| Barrel Outer Diameter  | 0.305 inches (7.75 mm)  |
|------------------------|-------------------------|
| Needle Inner Diameter  | 0.0059 inches (0.15 mm) |
| Needle Outer Diameter  | 0.0283 inches (0.72 mm) |
| Std Needle Dead Volume | 1.13 μL                 |
| Sales Quantity         | 1 syringe per pack      |

Specifications are subject to change without notice <a href="https://www.hamiltoncompany.com/laboratory-products/syringes/80930">https://www.hamiltoncompany.com/laboratory-products/syringes/80930</a>

**Spec. Version B** 



Visit <a href="https://www.hamiltoncompany.com/support">www.hamiltoncompany.com/support</a> to find a representative in your area or call us:

US: +1-775-525-3000 EU: +40-356-635-055